

## MEMO

DATE: March 14, 2020  
TO: Ryan Hammonds, Filbey Lab, Center for Brain Health  
FROM: Evan Kersey, Filbey Lab, Center for Brain Health  
SUBJECT: Freezer Alarm

### Introduction

There have been multiple issues related to the sample freezer recently where there has been a failure to notice or log when the temperature fluctuates outside safe ranges. A solution is a freezer alarm allowing immediate notification if the temperature fluctuates outside a certain safe range ( $T > -75^{\circ}\text{C}$ ). Building one from scratch will allow us to certify the probe will meet our extreme temperature requirements while also minimizing cost.

### Requirements

- I. Mitigate physical harm to the freezer
  - A. Mitigate any additional risk to samples in the freezer
  - B. Consult HVAC expert prior to any risky intervention
- II. Reliable
  - A. Handle power supply voltage fluctuation due to the compressor on the freezer kicking on
  - B. Provide notification if power is disrupted on the refrigerator circuit
  - C. All parts must be certified to work in a given temperature range
- III. Log temperature over time
  - A. Down time due to upgrades should be limited to 30 minutes per day
- IV. Alert lab personnel immediately if temperature falls outside of safe range
  - A. Notify personnel offsite
  - B. Notify personnel on site

### Proposed Solution

1. Insert a cryoprobe through the gasket into the freezer compartment.
  - a. This satisfies requirement I as it is a noninvasive solution. The diameter of the cable should be minimized.
  - b. This step must be successful in order to move forward with any other step.
  - c. Procedure.
    - i. Mount probe in freezer using glue or tape.
    - ii. Apply gentle heat to a 1in portion of the gasket to soften the rubber where the cable must pass through. This will help the gasket form a better seal.
    - iii. Close the freezer door and check for any leaks with the temp gun.

- d. If the procedure is successful, order other required parts.
- 2. Connect cryoprobe to an analog to digital converter chip.
  - a. The signal that comes through the wire from the probe is a very weak signal, too weak to be read by standard computer ports.
  - b. There are a wide range of price points here with different functionalities. If the selected component fails it is easily replaceable.
- 3. Connect converter to a control board.
  - a. There may be additional cables needed based on the exact model selected but they should be relatively inexpensive.
- 4. Write software logging temperature over time on a local thumb drive.
  - a. Proof of concept.
  - b. Allows verification that the components work with an acceptable degree of accuracy via the digital display on the freezer.
- 5. Interface the control board remotely.
  - a. Interfacing with a remote sensor could be difficult due to cometNet restrictions. Work with IT to allow direct access to the local server.(bhcent01)
  - b. Software on the server can pull data from the sensor and do something with it. Ideas include the following:
    - i. Chatbot on slack channel giving temperature updates every 30 minutes for peace of mind.
      - 1. Chatbot allows clients to query information about the sensor.
      - 2. I have a chatbot software written in Java (not the lab standard) and am able to convert it to pull information from a sensor instead of a RESTful API.

## Parts List

Cryoprobe/Thermocouple*	\$43.95 (-100C, 1200C) Type K [4]
Thermocouple -> male adapter	\$6.95[4]
Analog to Digital converter (compatible with thermocouple type)	\$14.95 [5][6]
Control Board, arduino uno rev3	\$22.00[7]
Surge protector	11.99[10]
glue/tape rated for -80C	
Temperature measurer tool (HVAC tech will have one if maintenance is done)	
Raspberry pi 4**	\$35[9]
Additional misc parts (breadboard cables, sd card, device casing)	<\$30
total	\$152.85

\*buy this first, if the cable is too thick to fit through the gasket then another approach is needed

\*\*this may be necessary to link to link to bhcent01 server more easily and also gives local ssh access

\*\*\*prices do not include shipping or tax

\*\*\*\*prices collected on 03/14/2020

## Alternative solutions

While alternative solutions do exist [8] most are unsuitable for our needs due to the probe and analog to digital converter not being rated to work in -80C conditions while also costing significantly more. This could be bypassed by buying an additional probe that is certified to work in these conditions however this would require doing a large amount of work on an undocumented commercial product. Costs associated with the existing alarms are around \$250 for hardware and \$250 for software a total of \$500 not including installation.[8] I was unable to find an option that met the -80C requirement out of the box.

## References

[1]<https://opensource.com/article/17/12/how-build-custom-iot-hardware-arduino>

Outlines procedure for creating an IOT sensor

[2]<https://www.mccdaq.com/TechTips/Measuring-Thermocouples-with-Raspberry-Pi.aspx>

Outlines requirement for hardware converting analog to digital board

The board for the RPI costs too much so going with an arduino and a cheap chip.

The arduino will need a breadboard

[https://www.youtube.com/watch?v=0KJ8H\\_SUWp0](https://www.youtube.com/watch?v=0KJ8H_SUWp0)

Video depicts arduino project with temp probe sensor

[4]<http://thesensorconnection.com/thermocouples/all-thermocouples/weld-pad-surface-temperature-thermocouple-probe>

thermocouple price

[5]<https://www.adafruit.com/product/269>

Digitizer price

[6]<https://cdn-shop.adafruit.com/datasheets/MAX31855.pdf>

Digitizer technical information

[7]<https://store.arduino.cc/usa/arduino-uno-rev3>

Arduino price

[8]<https://www.controlbyweb.com/applications/freezer-monitoring.html>

Alternative freezer alarm

[9]<https://www.raspberrypi.org/products/raspberry-pi-4-model-b/>

Rpi price, may become necessary to connect arduino to the bhcent server and add local functionality

[10][https://www.amazon.com/AmazonBasics-6-Outlet-Protector-2-Pack-2-Foot/dp/B014EKQ5AA?ref\\_=s9\\_apbd\\_otopr\\_hd\\_bw\\_b3C6a&pf\\_rd\\_r=E28ZA7MMR4GV74XJNHBH&pf\\_rd\\_p=c07a3cb9-6785-5a70-99d4-77ad257ef2c5&pf\\_rd\\_s=merchandised-search-10&pf\\_rd\\_t=BROWSE&pf\\_rd\\_i=761520](https://www.amazon.com/AmazonBasics-6-Outlet-Protector-2-Pack-2-Foot/dp/B014EKQ5AA?ref_=s9_apbd_otopr_hd_bw_b3C6a&pf_rd_r=E28ZA7MMR4GV74XJNHBH&pf_rd_p=c07a3cb9-6785-5a70-99d4-77ad257ef2c5&pf_rd_s=merchandised-search-10&pf_rd_t=BROWSE&pf_rd_i=761520)